

Cavitation/RF/Vacuum/Lipolaser



Principle of treatment

The Cavitation&Vacuum ---slimming, body and face-shaping:
Vacuum supersonic fat-explosion machine can promote tissue metabolism,repel the cellulite, which make fat granule in different depths and shocking make fusion energy, and the deepest receptor can get granule in 3CM thick of skin, and rapid expansion and the rapid heating of the cell wall that exceed the elastic limits and arises cells broken, so that the combustion of fatty acids at the same time broken out of a pot, not only the burning consumption of fatty acids but also more directly emptying so rapidly reduce the size of granule.



The RF---skin tightening, wrinkle removal:

High frequency wave permeates epidemis and take effect on the collagen-rich dermal tissue directly which makes water molecules in the skin move to and fro widely. The move produces controllable heat and results in the contraction of collagen and generation of new collagen and dermal remodeling. As a result, the skin looks tighter, wrinkles reduce gradually

The Lipo laser---weight loss

It emits low levels of laser energy, which creates a chemical signal in the fat cells, breaking down the stored triglycerides into free fatty acids and glycerol and releasing them though channels in the cell membranes. The fatty acids and glycerol are then transported around the body to the tissues that will use them during metabolism to create energy. This process of fatty acid release is a natural response of the body when the body needs to used stored energy reserves, thus Lipo is not creating any unnatural reaction in the body nor does it affect or damage any surrounding structures such as skin, blood vessels and peripheral nerves. A period of exercise post treatment will ensure the complete metabolism and thus elimination from the body of the freed fatty acids.

Application

- 1. Fat burning, Slimming, Body shaping
- 2. Effective promoting tissue metabolism and blood circulation, better for whiten skin
- 3. Improve orange peel organization
- 4. Strengthen the skin elasticity
- 5. Repair stretch marks
- 6. Anti-aging for face and body

Features

- 1. Adopt the leading-edge technology ultrasonic lipolysis in the world
- 2. Color touch screen with clear treatment parameters, friendly operation
- 3. Digital frequency control system, energy output evenly,precisely and steady, make sure parameter can't drifting
- 4. Suitable for all kinds of skin
- 5. 5--- 20 minutes per session treatment, no impact on your job or study
- 6. No pain, no cut, no scar, no any side-effect.
- 7. No rebound, long-lasting result!
- 8. Easy operation, easy study.

Treated Scope

- Vasodilator
- Fatness of orange skin
- Tighten the skin tissue, prevent slack
- Lipolysis
- Enhance Flexibility of the fibroblast
- Edema obesity
- Solid obesity
- Ecomposition of the fatness on the whole or partial body.

MACHINE INSTALLATION

		Lipo laser Wave lenght:635nm-650nm Max output:350mw Mix output:5mw Frequency:1hz-1Mhz
		Tripolar RF Supply output:66V Frequency:600khz Power:35W
		Body RF (Six pole RF) Supply output:66V Frequency:4MHZ Power:40W
		Face RF (Four pole RF) Supply output:66V Frequency:1MHZ Power:35W
		Vacuum RF Supply output:66V Frequency:1MHZ Power:95W Air pressure:□80kpa Put pressure:>250kpa Airflow:>10L/min Noise level:70dB(30cm outside)
		40K Cavitation Supply output:150V Frequency:40KHZ Power:25W